


# Mapping the test results from cobas b 101 to create automatic screening terms in Medtech.


When the test result comes into Medtech and it is filed, we need to make sure it creates an automatic screening term of the test result that came via HealthLink. **There will be an Hba1c test and a Lipids result which will have 5 tests within it that need to be mapped.**

Please note this will have to be completed prior to filing the test result within Medtech

1. Once you have logged into Medtech you need to click in on setup> Inbox/Outbox> Inbox screening.




2. Search for the PREDICT facility id of Hba1C and then double click into it.



Facility	Subject Description	Result Description	Screening Term	Screening Value
PREDICT	Hba1C	HBA1C		
PREDICT	Hdl	HDL		
PREDICT	Ldl	LDL		
PREDICT	Non-Hdl	NON-HDL		
PREDICT	Tc-Hdl	TC-HDL		
PREDICT	Tot-Chol	TOT-CHOL		
PREDICT	Tri	TRI		

3. Under 'Screening Term', double click and search for your Hba1C screening term and the same for 'Screening Value', you need to choose 'Hba1C' and press OK when finished.



**HBA1C (HBA1C)**

Main | **Audit**

Inbox Result

Facility: PREDICT

Subject Code: PREDICT

Subject Description: Hba1C

Result Code: HBA1C


Result Description: HBA1C

Screening Term

Screening Term: **HbA1c (HBA1C)**


Screening Value: Head Circumference (HCIR)  
Health Check (HC)  
Hearing (HEAR)  
Height (HE)  
Ideal body weight (IBW)

Inactive:




4. Look for the 'PREDICT' Facility ID with the subject description 'HDL' and the result description 'HDL' and double click into it. Under 'Screening Term' you need to choose your screening term that you use to save Lipid results.

*Note – Subject description may differ to the one shown below*



5. For 'Screening Value' you need to choose HDL Cholesterol and then press OK when finished.



**View Inbox Screening**

**HDL (HDL)**

Main | Audit

Inbox Result

Facility: PREDICT

Subject Code: PREDICT

Subject Description: Hdl

Result Code: HDL

Result Description: HDL

Screening Term


Screening Term: Lipids (LIPS)

Screening Value: [Dropdown menu open]

Inactive:

- Total Chol
- Triglycerides
- HDL Cholesterol**
- LDL Cholesterol

- Look for the 'PREDICT' Facility ID with the subject description 'LDL' and the result description 'LDL' and double click into it. Under 'Screening Term' you need to choose your screening term that you use to save Lipid results. For 'Screening Value' you need to choose 'LDL' Cholesterol'.



**View Inbox Screening**

**LDL (LDL)**

Main | Audit

Inbox Result

Facility: PREDICT

Subject Code: PREDICT

Subject Description: Ldl

Result Code: LDL

Result Description: LDL


Screening Term

Screening Term: Lipids (LIPS)

Screening Value: [Dropdown menu open]

Inactive:

- Lipids (LIPS)
- Liver Function Test (LFT)
- Mammography (MY)
- Maternity (MAT)
- Maternity 1 (MAT1)
- Maternity 2 (MAT2)



**View Inbox Screening**

**LDL (LDL)**

Main | Audit

Inbox Result

Facility: PREDICT

Subject Code: PREDICT

Subject Description: Ldl

Result Code: LDL

Result Description: LDL

Screening Term


Screening Term: Lipids (LIPS)

Screening Value: LDL Cholesterol

Inactive:

- Total Chol
- Triglycerides
- HDL Cholesterol
- LDL Cholesterol
- TC/HDL Ratio

7. Look for the 'PREDICT' Facility ID with the subject description 'TC-HDL' and the result description 'TC-HDL' and double click into it. Under 'Screening Term' you need to choose your screening term that you use to save Lipid results. For 'Screening Value' you need to choose 'TC/HDL Ratio'.



**View Inbox Screening**

**TC-HDL (TC-HDL)**

Main | Audit

Inbox Result

Facility: PREDICT

Subject Code: PREDICT

Subject Description: Tc-Hdl

Result Code: TC-HDL

Result Description: TC-HDL


Screening Term

Screening Term: Lipids (LIPS)

Screening Value: Lipids (LIPS)

Inactive:

- Liver Function Test (LFT)
- Mammography (MY)
- Maternity (MAT)
- Maternity 1 (MAT1)
- Maternity 2 (MAT2)



**View Inbox Screening**

**TC-HDL (TC-HDL)**

Main | Audit

Inbox Result

Facility: PREDICT

Subject Code: PREDICT

Subject Description: Tc-Hdl

Result Code: TC-HDL

Result Description: TC-HDL

Screening Term


Screening Term: Lipids (LIPS)

Screening Value:

Inactive:

- Total Chol
- Triglycerides
- HDL Cholesterol
- LDL Cholesterol
- TC/HDL Ratio**

- Look for the 'PREDICT' Facility ID with the subject description 'TOT-CHOL' and the result description 'TOT-CHOL' and double click into it. Under 'Screening Term' you need to choose your screening term that you use to save Lipid results. For 'Screening Value' you need to choose 'Total Cholesterol'.



**View Inbox Screening**

**TC-HDL (TC-HDL)**

Main | Audit

Inbox Result

Facility: PREDICT

Subject Code: PREDICT

Subject Description: Tc-Hdl

Result Code: TC-HDL

Result Description: TC-HDL


Screening Term

Screening Term: Lipids (LIPS)

Screening Value: Lipids (LIPS)

Inactive:

- Liver Function Test (LFT)
- Mammography (MY)
- Maternity (MAT)
- Maternity 1 (MAT1)
- Maternity 2 (MAT2)



**View Inbox Screening**

**TOT-CHOL (TOT-CHOL)**

Main | Audit

Inbox Result

Facility: PREDICT

Subject Code: PREDICT

Subject Description: Tot-Chol

Result Code: TOT-CHOL

Result Description: TOT-CHOL

Screening Term


Screening Term: Lipids (LIPS)

Screening Value: [Dropdown]

Inactive:

- Total Chol
- Triglycerides
- HDL Cholesterol
- LDL Cholesterol

- Look for the 'PREDICT' Facility ID with the subject description 'TRI' and the result description 'TRI' and double click into it. Under 'Screening Term' you need to choose your screening term that you use to save Lipid results. For 'Screening Value' you need to choose 'Triglycerides'. Press OK



**View Inbox Screening**

**TRI (TRI)**

Main | Audit

Inbox Result

Facility: PREDICT

Subject Code: PREDICT

Subject Description: Tri

Result Code: TRI

Result Description: TRI

Screening Term

Screening Term: Lipids (LIPS)

Screening Value: Lipids (LIPS)

Inactive:

- Lipids (LIPS)
- Liver Function Test (LFT)
- Mammography (MY)
- Maternity (MAT)
- Maternity 1 (MAT1)
- Maternity 2 (MAT2)

**View Inbox Screening**

**TRI (TRI)**

Main | Audit

Inbox Result

Facility: PREDICT

Subject Code: PREDICT

Subject Description: Tri

Result Code: TRI

Result Description: TRI

Screening Term

Screening Term: Lipids (LIPS)

Screening Value: Triglycerides

Inactive:  Triglycerides

HDL Cholesterol

LDL Cholesterol

The mapping is now complete.